



Blood Collection Set

 Building a Healthier
Tomorrow

RAMY

1st Industrial City | P.O. Box 7917 |
Road 22|Cross 11 | Dammam-31472|
Kingdom of Saudi Arabia

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al-Shifa supplies Blood Collection Set also called scalp vein set or winged infusion set, is a medical device for either withdrawing blood from a vein of patient for diagnosis or administering IV therapy into a vein.

The design of Blood Collection Set consists of an ultra-sharp siliconized needle to minimize penetration and gliding force as well as patient discomfort. The thin wall of the needle provides better flow rate per gauge size. The wings of the Blood Collection set are flexible and easy to grip and attachment with the skin.

The wings of the Blood Collection set are color coded for instant identification of needle size.

Al-Shifa Medical Co. supplies different types and sizes of Blood Collection Set.

Blood Collection Set

Product Description

Sterile, Single Use – Blood Collection Set consists of a hub, needle, cannula of penetrating, tube with rubber sleeve and a needle protector, where necessary venous accessories. The needle is used matching with vacuum blood collection tube mainly used for collection human venous blood or plasma. When using it, the operator of the device should be a trained professional medical staff.

Product Intended Use

Sterile, Single Use – Blood Collection/Infusion Set is used for blood collection and/or the short-term infusion of intravenous fluids. The winged needle is designed with a safety shield, which can be activated to cover the needle immediately following venipuncture to aid in the protection against accidental needlestick injury.

Product Features

- Sharp Stainless-Steel Needle for atraumatic insertion
- Thin walled for optimal flowrate preserving blood cell integrity
- Ergonomic grip, one handed activation of safety mechanism
- Colour Coded for easy needle size identification
- Free from Air and Liquid leakage
- Pyrogen Free
- Single Use, Disposable
- Sterile
- Use by Healthcare Professionals only

Product Types and Variants

Blood Collection Set with Pre-Attached Holder



Size	Color Code	Catalogue #
G: 19 x 3/4" x 7"		F35111975070
G: 21 x 3/4" x 7"		F35112175070
G: 23 x 3/4" x 7"		F35112375070
G: 25 x 3/4" x 7"		F35112575070
G: 27 x 3/4" x 7"		F35112775070

Blood Collection Set (Regular Scalp vein)



Blood Collection Set with Plug Cap (Safety Scalp Vein)



Size	Colour Code	Catalogue #
G: 18 x 3/4"		F05041819300
G: 19 x 3/4"		F05041919300
G: 20 x 3/4"		F05042019300
G: 21 x 3/4"		F05042119300
G: 22 x 3/4"		F05042219300
G: 23 x 3/4"		F05042319300
G: 25 x 3/4"		F05042519300
G: 25 x 3/4"		F05042519200
G: 27 x 3/4"		F05042719300

Size	Colour Code	Catalogue #
G: 18 x 3/4"		F05031819300
G: 19 x 3/4"		F05031919300
G: 20 x 3/4"		F05032019300
G: 21 x 3/4"		F05032119300
G: 22 x 3/4"		F05032219300
G: 23 x 3/4"		F05032319300
G: 24 x 3/4"		F05032419300
G: 25 x 3/4"		F05032519300
G: 26 x 3/4"		F05032619300
G: 27 x 3/4"		F05032719300

Blood Collection Set with Luer Adaptor and BCN (Safety) w/o holder



Size	Color Code	Catalogue #
G: 18 x 3/4"		F33351819191
G: 19 x 3/4"		F33351919191
G: 20 x 3/4"		F33352019191
G: 22 x 3/4"		F33352219191
G: 23 x 3/4"		F33352319191
G: 25 x 3/4"		F33352519191
G: 27 x 3/4"		F33352719191

Blood Collection Set with Luer Adaptor and BCN (Regular) w/o holder



Size	Color Code	Catalogue #
G: 18 x 3/4"		F33351819190
G: 19 x 3/4"		F33351919190
G: 20 x 3/4"		F33352019190
G: 21 x 3/4"		F33352119190
G: 22 x 3/4"		F33352219190
G: 23 x 3/4"		F33352319190
G: 25 x 3/4"		F33352519190
G: 27 x 3/4"		F33352719190

Multi-Sample Blood Collection Needle



Size	Catalogue #
G: 20 x 1"	F80012010000
G: 20 x 1 1/2"	F80012015000
G: 21 x 1"	F80012110000
G: 21 x 1 1/2"	F80012115000
G: 22 x 1"	F80012210000
G: 22 x 1 1/2"	F80012215000
G: 23 x 1"	F80012310000
G: 23 x 1 1/2"	F80012315000

Our Company



Our passion is excellence, reliability, and creating innovative solutions. Since our establishment in 1980, we have committed to ourselves to making the healthcare environment safe and clean by manufacturing the latest and highest-quality medical disposable. We achieved that through a dedicated and integrated team that strives to make our company a leader in the field. We have achieved many major accomplishments in the past, including owning the largest medical injection syringe factory in the Gulf and contributing to the fight against the coronavirus pandemic, by manufacturing over 60 million vaccine doses for Saudi Arabia and Dubai. With many other promising projects and accomplishments that have made us a focus of attention and a hub for creative partners, we aim to be the first entity in our field, and we spare no effort in continuous development and innovation.