

Insulin Pen Needle

 **Building a Healthier
Tomorrow**

Inspire

1st Industrial City | P.O. Box 7917 |
Road 22 | Cross 11 | Dammam-31472 |
Kingdom of Saudi Arabia

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Insulin Pen Needles are used for the treatment of Diabetes Patients.

Al Shifa Medical Co. Supplies Insulin Pen Needles compatible with the pen injectors such as NovoPen®5, Autopen®24, Gansulin Pen, Ypsopen®, UNIPEN®, XiuLinPen, HumaPen® Ergo II, Byetta, Victoza Pen.

The needle of the insulin pen is designed considering thin wall technology. The thin wall of the needle supports comfortable skin penetration which minimizes the patient's discomfort during injections.

Insulin Pen Needle

Product Description

Sterile, Single Use – Insulin Pen Needles are attached with compatible Insulin Pen Injectors which is used for subcutaneous injection of insulin in the treatment of diabetes. Triple bevel needle tip and silicone lubrication ensure pain Free injection. The thin wall technology of needles improves the drug flow resulting in less thumb force and time to deliver the dose.

Product Intended Use

Sterile, Single Use – Insulin Pen Needle is intended for use with Insulin pen injector devices for the subcutaneous injection of insulin in diabetic patient.

Product Features

- Compatible with Branded Insulin Pen Injectors
- Sharp, Stainless Steel thin wall needle technology
- Design provides strong gripping during connection with Insulin Pen Injectors.
- Pyrogen Free
- Single Use, Disposable
- Sterile
- Use by Healthcare Professionals only


Product Types and Variants

Insulin Pen Needle

Size	Catalogue #
G: 31 x 6mm	F05103106000
G: 31 x 5mm	F05103105000
G: 31 x 4mm	F05103104000
G: 32 x 6mm	F05103206000
G: 32 x 5mm	F05103205000
G: 32 x 4mm	F05103204000
G: 33 x 6mm	F05103306000
G: 33 x 5mm	F05103305000
G: 33 x 4mm	F05103304000



Our Company



Our passion is excellence, reliability, and creating innovative solutions.

Since our establishment in 1980, we have committed to ourselves to making the healthcare environment safe and clean by manufacturing the latest and highest-quality medical disposable.

We achieved that through a dedicated and integrated team that strives to make our company a leader in the field.

We have achieved many major accomplishments in the past, including owning the largest medical injection syringe factory in the Gulf and contributing to the fight against the coronavirus pandemic, by manufacturing over 60 million vaccine doses for Saudi Arabia and Dubai.

With many other promising projects and accomplishments that have made us a focus of attention and a hub for creative partners, we aim to be the first entity in our field, and we spare no effort in continuous development and innovation.

Inspire

