



I.V Cannula with Safety

With & Without
Injection port

 **Building a Healthier
Tomorrow**

RAMY

1st Industrial City | P.O. Box 7917 |
Road 22 | Cross 11 | Dammam-31472 |
Kingdom of Saudi Arabia

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Safety I.V Cannula is one of the most common device use in infusion therapy.

Al Shifa Medical Co. Supplies Safety I.V Cannula which is designed for accuracy and comfort in infusion therapy of patients as well as the product provides a Safety Feature which secure the patients and Healthcare professionals from needle stick injuries.

Safety I.V Cannula have an ULTRA-SHARP tip which minimizes patient discomfort during Infusion Therapy along with its High-Quality kink resistant FEP radio-opaque catheter for easy insertion with optimal flowrate.

The Safety I.V Cannula's are with color coded Hub which indicates its Gauge size.

Safety I.V Cannula

with & without Injection Port

Product Description

Intra Venous Cannula with safety feature is used for intravenous infusion of parenteral fluid /Medication/ withdrawal of blood sample take place in a patient. It comprises of several parts, body in which catheter is mounted, a needle which is mounted on a hub, and it passes through the body and catheter. It is used as a introduce needle catheter glide and enter the vein. There is a hub cover which is mounted on needle hub, it is designed in such a way that air is displaced when the needle tip enters the blood path the only part which remain in contact with patient is catheter, catheter mount and wing body.

Product Intended Use

Sterile Single Use – A peripheral intravenous catheter (PIVC) with safety feature is a thin plastic tube inserted into a vein using a needle. PIVCs allow for the administration of medications, fluids and/or blood products.

Product Features

- Sharp and siliconized with safety feature stainless-steel needle for smooth insertion
- Clear and Transparent Flexible TEFLON Catheter
- With/without Flexible wings
- With/without Injection Flip Port Cap
- With Protective Needle Cover to and Safety feature prevent sharp injuries
- Accurate in Dimensions
- Colour coded catheter hub and flip port cap for easy needle size identification.
- Free from Air and Liquid leakage
- Radio-Opaque
- Pyrogen Free
- Single Use, Disposable
- Sterile
- For Healthcare Professionals use only

Product Types and Variants

Safety I.V Cannula with Wings with Injection Port



Size	Color Code	Catalogue #
G: 14 x 45mm		F22210121400
G: 16 x 45mm		F22210121600
G: 17 x 45mm		F22210121700
G: 18 x 45mm		F22210121800
G: 18 x 30mm		F22210121800
G: 20 x 32mm		F22210122000
G: 22 x 25mm		F22210122200
G: 24 x 19mm		F22210122400
G: 26 x 19mm		F22210122600

Safety I.V Cannula with Wings & without Injection Port- Blood Control



Size	Color Code	Catalogue #
G: 14 x 45mm		F22210001445
G: 16 x 45mm		F22210001645
G: 16 x 30mm		F22210001630
G: 18 x 30mm		F22210001830
G: 18 x 45mm		F22210001845
G: 20 x 30mm		F22210002030
G: 22 x 25mm		F22210002225
G: 24 x 19mm		F22210002419

Safety I.V Cannula

with & without Injection Port

Safety I.V Cannula without Wings & without
Injection Port and without blood collection
– Pen Type



Size	Colour Code	Catalogue #
G: 14 x 30mm		F22210131430
G: 14 x 45mm		F22210131445
G: 16 x 45mm		F22210131645
G: 18 x 45 mm		F22210131845
G: 20 x 25mm		F22210132025
G: 20 x 32mm		F22210132032
G: 22 x 25mm		F22210132225
G: 24 x 19mm		F22210132419
G: 26 x 19mm		F22210132619

Safety I.V Cannula without Wings & without
Injection Port- Blood Control



Size	Colour Code	Catalogue #
G: 14 x 45mm		F22210141445
G: 16 x 45mm		F22210141645
G: 16 x 30mm		F22210141630
G: 18 x 30mm		F22210141830
G: 18 x 45mm		F22210141845
G: 20 x 30mm		F22210142030
G: 22 x 25mm		F22210142225
G: 24 x 19mm		F22210142419
G: 26 x 19mm		F22210102619

Safety IVC with wings & without injection port with
blood control/notched needle, push button



Size	Colour Code	Catalogue #
G: 14 x 30mm		F22210001430
G: 14 x 45mm		F22210001445
G: 16 x 30mm		F22210001630
G: 16 x 45 mm		F22210001645
G: 18 x 30mm		F22210001830
G: 18 x 32mm		F22210001832
G: 18 x 45mm		F22210001845
G: 20 x 32mm		F22210002032
G: 20 x 30mm		F22210002030
G: 22 x 25mm		F22210002225
G: 24 x 19mm		F22210002419
G: 26 x 19mm		F22210002619

RAMY