



RAMY

Hypodermic Syringe

(With Needle)

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al Shifa Medical Co. Supplies Hypodermic Syringes which are designed for accuracy and comfort in treatments of patients.

Hypodermic Syringe with needle produce both Luer Lock and Luer Slip designs. Its Clear Barrels with bold scale markings permit quick visual inspection for precise dosage.

An integrated plunger stop helps prevent the plunger from being accidentally pulled out and its textured thumb press design reduce the risk of slippage during administration.

The shape and size of flange provide stability and comfort during use.

1st Industrial City | P.O. Box 7917 |
Road 22 | Cross 11 | Dammam-31472 |
Kingdom of Saudi Arabia

Product Description

Sterile Single Use –Ramy Hypodermic Syringe with needle consist of Core Plunger, Syringe Barrel, Plunger Rod, Stainless Steel Needle and Protective Cap.

Product Intended Use

Sterile Single use Hypodermic Syringe with Needle is intended to be used for medical purposes to inject fluid into or withdraw fluid from body.

Product Features

- Clear and Transparent Barrel for Easy Visualization of medication/fluids.
- Fixed printed Graduation Scale
- Smooth movement of plunger requiring minimum injecting force and withdrawal force.
- Free from liquid and Air leakage.
- Sharp needle and smooth surface with siliconization for easy penetration.
- Sterile by E.O.
- Pyrogen Free.
- Single Use

Product Types and Variants

Hypodermic Syringe Luer Lock with Needle

Size	Needle Gauge	Catalogue #
3ml	23G x 1"	F 30302310220
	23G x 1-¼"	F 30302312220
5ml	23G x 1"	F 30502310220
	23G x 1-½"	F30502315220
	22GX1-¼"	F 30502212220
10ml	21GX1-½"	F 30502115220
	23G x 1"	F 30102310220
	22GX1-¼"	F 30102212220
20ml	21GX1-½"	F 30102115220
	23G x 1"	F 20002115220
30ml	21GX1-½"	F 20002115220
50ml	21GX1-½"	F 50002115220



Hypodermic Syringe Luer Slip with Needle

Size	Needle Gauge	Type	Catalogue #
3ml	23G x 1-¼"	CENTRIC	F 30302312330
	23G x 1"	CENTRIC	F 30302310330
5ml	23GX1"	CENTRIC	F 30502310330
	22GX1-¼"	ECCENTRIC	F 31102212330
10ml	21GX1-½"	ECCENTRIC	F 31102115330
	21G x 1-1/2"	CENTRIC	F 30102115330
20ml	21GX1-½"	ECCENTRIC	F 20002115330



RAMY Hypodermic Syringes (With Needle)

RAMY Hypodermic Syringes (With Needle)

Vaccination Syringe with Fixed Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F 13002606440
	25G x 5/8"	F 13002506440
	25Gx1	F 13002510 440

Vaccination Syringe Luer Lock with Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F 13002606220

Tuberculin Syringe Luer Slip with Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F 13012606330
	25G x 5/8"	F 13012506330

Vaccination Syringe Luer Slip with Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F 13002606330
	25G x 5/8"	F 13002506330

Tuberculin Syringe with Fixed Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F 13012606440
	25G x 5/8"	F 13012506440

Tuberculin Syringe with Luer lock Needle

Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F1301260633

