

RAMY Hypodermic Syringe (Without Needle)

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al Shifa Medical Co. Supplies Hypodermic Syringes which are designed for accuracy and comfort in treatments of patients.

Hypodermic Syringe without needle produce both Luer Lock and Luer Slip designs. Its Clear Barrels with bold scale markings permit quick visual inspection for precise dosage.

An integrated plunger stop helps prevent the plunger from being accidentally pulled out and its textured thumb press design reduce the risk of slippage during administration.

The shape and size of flange provide stability and comfort during use.

Product Description

Sterile Single Use – Ramy Hypodermic Syringe with needle consist of Core Plunger, Syringe Barrel, and Plunger Rod.

Product Intended Use

Sterile Single use – Hypodermic Syringe without Needle is intended to be used for medical purposes to inject fluid into or withdraw fluid from body.


Product Features

- Clear and Transparent Barrel for Easy Visualization of medication/fluids.
- Fixed printed Graduation Scale
- Smooth movement of plunger requiring minimum injecting and/or withdrawal force.
- Free from liquid and Air leakage.
- Sterile by E.O.
- Pyrogen Free
- Latex Free
- Single Use

Product Types and Variants

Hypodermic Syringe – Luer Lock

Size	Catalogue #
3 ml	F30300000220
5 ml	F30500000220
10 ml	F30100000220
20 ml	F20000000220
30 ml	F30000000220
50 ml	F50000000220
60 ml	F60000000220



Hypodermic Syringe – Luer Slip

Size	Type	Catalogue #
2.5 ml	-	F30252115330
3 ml	Centric	F30300000330
5 ml	Centric	F30500000330
10ml	Centric	F31100000330
10 ml	Centric	F30100000330
20 ml	Eccentric	F20000000330
50 ml	ECCENTRIC / SIDE NOZZLE)	F50000000330
50 ml	Catheter Tip	F50000000660



Hypodermic Syringes (Without Needle)