



Urine Drainage Bag

 **Building a Healthier
Tomorrow**
OInspire

1st Industrial City | P.O. Box 7917 |
Road 22 | Cross 11 | Dammam-31472 |
Kingdom of Saudi Arabia

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al-Shifa Medical Co. supplies urine collection bags which are designed to collect urine drained from the bladder via a catheter or sheath.

Al-Shifa Medical Co. provides Normal Urine Bag, Urine drainage bag with sampling port, Urine bag for neonates and urine leg bags, these types of urine drainage bags are available in various sizes and styles to meet different needs.

Urine Leg Bags are usually made of waterproof materials and come with adjustable straps to ensure a secure connection.

Urine Bags are typically hung on a stand or bed frame and come with a secure closure system to prevent overflow and leakage.

Urine Bags are an essential component of catheter care. They provide a convenient and hygienic way to collect and store urine, which helps to reduce infection risk and improve patient comfort.

Urine Drainage Bag

Product Description

Sterile, Single Use – Urine Drainage Bag consist of leakage free bag, soft and kink resistant tube which provides route to transfer the urine and connectors to make feasible secure connections between urine bag and catheter.

Product Intended Use

Sterile, Single Use – Urine Drainage Bag is intended to connect with indwelling catheter for urine collection and drainage.

Product Features

- With large volume capacity
- With different lengths of tubing with bonded connector
- Non-return valve
- Integrated Hanger with rope for easy attachment to the desired location
- Clear printed incremental and accurate measuring marking
- Non-Toxic, Pyrogen Free and Latex-Free
- Single Use, Disposable
- Sterile
- Use by Healthcare Professionals only

Product Types and Variants

Urine Drainage Bag – Regular

Size	Tube Length	Catalogue #
3000ml	50"	FUB130005006
2800ml	50"	FUB128005006
2500ml	50"	FUB125005006
2200ml	50"	FUB122005006
2000ml	50"	FUB120005006
1500ml	42"	FUB115004211
1000ml	35"	FUB110003510
2000ml	90"	FUB200000120



Urine Drainage Bag – Regular

Size	Tube Length	Catalogue #
500ml	25"	FUB305002510
800ml	25"	FUB308002510
1000ml	25"	FUB310002510



Urine Drainage Bag – with Sampling Port



Size	Tube Length	Catalogue #
3000ml	50"	FUB230005006
2800ml	50"	FUB228005006
2500ml	50"	FUB225005006
2200ml	50"	FUB222005006
2000ml	50"	FUB220005006
1500ml	42"	FUB215004211
1000ml	35"	FUB210003510

Urine Drainage Bag - Pediatric



Size	Catalogue #
50ml	FUB00050BA
200ml	FUB2000500FC
100ml	FUB20001000H

Urine Drainage Bag – Fine chamber



Size	Catalogue #
50ml	FUB2000500URO

Luxury Urine Drainage Bag – Double Hook



Size	Tube Length	Catalogue #
2000ml	100cm	FUB200100HH

Urine Drainage Leg bag with Valve



Size	Tube Length	Catalogue #
500ml	25cm	FUB500100LEG
800ml	25cm	FUB800100LEG

Our Company



Our passion is excellence, reliability, and creating innovative solutions. Since our establishment in 1980, we have committed to ourselves to making the healthcare environment safe and clean by manufacturing the latest and highest-quality medical disposable. We achieved that through a dedicated and integrated team that strives to make our company a leader in the field. We have achieved many major accomplishments in the past, including owning the largest medical injection syringe factory in the Gulf and contributing to the fight against the coronavirus pandemic, by manufacturing over 60 million vaccine doses for Saudi Arabia and Dubai. With many other promising projects and accomplishments that have made us a focus of attention and a hub for creative partners, we aim to be the first entity in our field, and we spare no effort in continuous development and innovation.