



Olnspire Hypodermic Syringe

(With & Without
Needle)

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al Shifa Medical Co. Supplies Hypodermic Syringes which are designed for accuracy and comfort in treatments of patients.

Hypodermic Syringes with and without needle are produced in both Luer Lock and Luer Slip designs. Its Clear Barrels with bold scale markings permit quick visual inspection for precise dosage.

An integrated plunger stop helps prevent the plunger from being accidentally pulled out and its textured thumb press design reduce the risk of slippage during administration.

The shape and size of flange provide stability and comfort during use.

Product Description

Sterile Single Use – Inspire Hypodermic Syringe with & without needle consist of Core Plunger, Syringe Barrel, Plunger Rod, Stainless Steel Needle and Protective Cap.

Product Intended Use

Sterile Single use – Hypodermic Syringe with & without Needle is intended to be used for medical purposes to inject fluid into or withdraw fluid from body.

Product Features

- Clear and Transparent Barrel for Easy Visualization of medication/fluids.
- Fixed printed Graduation Scale
- Smooth movement of plunger requiring minimum injecting force and withdrawal force.
- Free from liquid and Air leakage.
- Sharp needle and smooth surface with siliconization for easy penetration.
- Sterile by E.O.
- Pyrogen Free.
- Single Use

Product Types and Variants

Hypodermic Syringe Luer Lock with Needle

Size	Needle Gauge	Catalogue #
3ml	22GX1-½"	F30312215330
	23G x 1-¼"	F30312312220
5ml	23G x 1"	F30512310220
	22GX1-¼"	F30512212220
10ml	21GX1-½"	F30512115220
	23G x 1"	F30112310220
20ml	22GX1-¼"	F30112212220
	21GX1-½"	F30112115220
30ml	21GX1-½"	F20012115220



Hypodermic Syringe Luer Slip with Needle



Size	Type	Needle Gauge	Catalogue #
3ml	Centric	23G x 1-¼"	F30312312330
		23G x 1"	F30312310330
5ml	-	23GX1.5"	F30302115220
	Centric	23GX1"	F30512310330
10ml	-	22GX1.5"	F30512315330
	Eccentric	22GX1-¼"	F31112212330
20ml	Eccentric	21GX1-½"	F31112115330
	Centric	21GX1-½"	F30112115330
	Eccentric	21GX1-½"	F20012115330

Hypodermic Syringes (With Needle)

Hypodermic Syringes (With Needle)

Tuberculin Syringe with Needle



Size	Needle Gauge	Type	Catalogue #
1ml	26G x 5/8"	Luer lock	F13022606220
	25G x 5/8"	Luer Slip	F13032506330
	25GX5/8"		F13022606220

Vaccination Syringe with Fixed Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F13022606440
	25G x 5/8"	F13022506440
	25G x 1"	F13032510440

Vaccination Syringe Luer Slip with Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F13022606330
	25G x 5/8"	F1302 2506 330

Vaccination Syringe Luer Lock with Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F13022606220

Tuberculin Syringe with Fixed Needle



Size	Needle Gauge	Catalogue #
1ml	26G x 5/8"	F13032606440

Hypodermic Syringes (With Needle)

Arterial Blood Sampling Syringe

Size	Needle Gauge	Preset Vol.	Heparin Lithium IU	Catalogue #
1ml	0.5 x 16mm	0.6ml	30IU	ARB1516
	0.5 x 25mm	0.6ml	30IU	ARB1525
	0.5 x 32mm	0.6ml	30IU	ARB1532
	0.6 x 16mm	0.6ml	30IU	ARB1616
	0.6 x 25mm	0.6ml	30IU	ARB1625
	0.6 x 32mm	0.6ml	30IU	ARB1632
	0.7 x 25mm	0.6ml	30IU	ARB1725
3ml	0.7 x 32mm	0.6ml	30IU	ARB1572
	0.5 x 16mm	1.6ml	80IU	ARB3516
	0.6 x 16mm	1.6ml	80IU	ARB3616
	0.6 x 25mm	1.6ml	80IU	ARB3625
	0.6 x 32mm	1.6ml	80IU	ARB3632
	0.7 x 25mm	1.6ml	80IU	ARB3725
	0.7 x 32mm	1.6ml	80IU	ARB3732




Arterial Blood Sampling Syringe – Safety



Size	Needle Gauge	Preset Vol.	Heparin Lithium IU	Catalogue #
1ml	0.5 x 16mm	0.6ml	30IU	ARB1516S
	0.5 x 25mm	0.6ml	30IU	ARB1525S
	0.5 x 32mm	0.6ml	30IU	ARB1532S
	0.6 x 16mm	0.6ml	30IU	ARB1616S
	0.6 x 25mm	0.6ml	30IU	ARB1625S
	0.6 x 32mm	0.6ml	30IU	ARB1632S
	0.7 x 25mm	0.6ml	30IU	ARB1725S
3ml	0.7 x 32mm	0.6ml	30IU	ARB1732S
	0.5 x 16mm	1.6ml	80IU	ARB3516S
	0.6 x 16mm	1.6ml	80IU	ARB3616S
	0.6 x 25mm	1.6ml	80IU	ARB3625S
	0.6 x 32mm	1.6ml	80IU	ARB3632S
	0.7 x 25mm	1.6ml	80IU	ARB3725S
	0.7 x 32mm	1.6ml	80IU	ARB3732S

Hypodermic Syringes

(Without Needle)



Hypodermic Syringe – Luer Lock Without Needle

Size	Catalogue #
3 ml	F30310000220
5 ml	F30510000220
10 ml	F30110000220
20 ml	F20010000220
30 ml	F30310000220
50 ml	F30310000220
60 ml	F60010000220

Hypodermic Syringe – Luer Slip Without Needle

Size	Type	Catalogue #
2.5 ml	-	F30262115330
3 ml	Centric	F30310000330
5 ml	Centric	F30510000330
10ml	Eccentric	F31110000330
	Centric	F30110000330
20 ml	Eccentric	F20010000330
50 ml	Eccentric	F50010000330
50 ml	Catheter Tip	F50010000660



Arterial Blood Sampling Syringe - Without Needle



Size	Preset Vol.	Heparin Lithium IU	Catalogue #
1ml	0.6ml	30IU	ARB1006
3ml	1.6ml	80IU	ARB3016

Vaccine Luer Lock – Without needle

Size	Catalogue #
1 ml	F13010000220

Tuberculin Luer slip – Without needle

Size	Catalogue #
1 ml	F13032606330

Our Company



Our passion is excellence, reliability, and creating innovative solutions.

Since our establishment in 1980's, we have committed to ourselves to making the healthcare environment safe and clean by manufacturing the latest and highest-quality medical disposable.

We achieved that through a dedicated and integrated team that strives to make our company a leader in the field.

We have achieved many major accomplishments in the past, including owning the largest medical injection syringe factory in the Gulf and contributing to the fight against the coronavirus pandemic, by manufacturing over 60 million vaccine doses for Saudi Arabia and Dubai.

With many other promising projects and accomplishments that have made us a focus of attention and a hub for creative partners, we aim to be the first entity in our field, and we spare no effort in continuous development and innovation.