



RAMY Hypodermic Needles

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al Shifa Medical Co. Supplies Hypodermic Needles for the Intramuscular, Intravenous and/or Intradermal for injecting of medicinal fluid or withdrawal of blood for diagnosis of patients.

Hypodermic Needles have an ULTRA-SHARP tip which minimizes patient discomfort during injection and/or withdrawal along with syringe and its High-Quality break resistance

The Hypodermic Needles are with color coded Hub which indicates its Gauge size.

Al-Shifa Medical also supply Hypodermic Safety Needle which minimizes the risk of sharp injuries to the medical staff as well as to the patients. Its innovative design enables the safety system with one click.

1st Industrial City | P.O. Box 7917 |
Road 22|Cross 11 | Dammam-31472|
Kingdom of Saudi Arabia

Product Description

Sterile Single Use – Hypodermic needles are normally made from a stainless-steel tube through a process known as tube drawing where the tube is drawn through progressively smaller dies to make the needle. The end is bevelled to create a sharp pointed tip letting the needle easily penetrate the skin.

Product Intended Use

Sterile Single Use – Hypodermic needle is used with syringes to inject substances or extract fluids from human body.

Product Features

- Sharpened tip for smooth.
- Siliconized for smooth penetration.
- Colour coded hub for easy Needle size identification.
- Latex and DEHP Free.
- Non-Toxic
- Pyrogen Free.
- Single Use
- Sterile by E.O.
- For Healthcare Professionals use only.

Product Types and Variants

Hypodermic Needle

Size	Color Code	Catalogue #
G: 16 x 1.5"		F05001615000
G: 18 x 1.5"		F05001815000
G: 19 x 1.25"		F05001912000
G: 19 x 1"		F05001910000
G: 19 x 1.5"		F05001915000
G: 20 x 1.5"		F05002015000
G: 21 x 1"		F05002110000
G: 21 x 1.5"		F05002115000
G: 22 x 1.5"		F05002215000
G: 22 x 1"		F05002210000
G: 22 x 1.25"		F05002212000
G: 23 x 1"		F05002310000
G: 23 x 1.5"		F05002315000
G: 23 x 1.25"		F05002312000
G: 24 x 1"		F05002410000
G: 25 x 1"		F05002510000
G: 25 x 1.5"		F05002515000
G: 25 x 5/8"		F05002506000
G: 26 x 5/8"		F05002606000
G: 26 x 0.5"		F05002605000
G: 27 x 0.5"		F05002705000
G: 30 x 0.5"		F05003005000



RAMY

Hypodermic Needles



Our Company



Our passion is excellence, reliability, and creating innovative solutions.

Since our establishment in 1980, we have committed to ourselves to making the healthcare environment safe and clean by manufacturing the latest and highest-quality medical disposable.

We achieved that through a dedicated and integrated team that strives to make our company a leader in the field.

We have achieved many major accomplishments in the past, including owning the largest medical injection syringe factory in the Gulf and contributing to the fight against the coronavirus pandemic, by manufacturing over 60 million vaccine doses for Saudi Arabia and Dubai.

With many other promising projects and accomplishments that have made us a focus of attention and a hub for creative partners, we aim to be the first entity in our field, and we spare no effort in continuous development and innovation.

Hypodermic Needle – with SAFETY

Size	Color Code	Catalogue #
G: 18 x 2"		F05011820000
G: 18 x 1 1/2"		F05011815000
G: 18 x 1 1/2"		F0501181500
G: 18 x 1"		F05011810000
G: 18 x 1"		F0501181000
G: 19 x 2"		F05011920000
G: 19 x 1 1/2"		F05011915000
G: 19 x 1"		F05011910000
G: 19 x 1"		F0501191000
G: 20 x 1 1/2"		F05012015000
G: 20 x 1 1/2"		F0501201500
G: 20 x 1"		F05012010000
G: 20 x 2"		F05012020000
G: 21 x 2"		F05012120000
G: 21 x 1 1/5"		F05012115000
G: 21 x 1 1/5"		F0501211500
G: 21 x 1"		F05012110000
G: 22 x 1 1/5"		F05012215000
G: 22 x 1 1/5"		F0501221500
G: 22 x 1 1/4"		F05012212000
G: 22 x 1"		F05012210000
G: 23 x 1 1/5"		F05012315000
G: 23 x 1 1/5"		F0501231500
G: 23 x 1 1/4"		F05012312000
G: 23 x 1"		F05012310000
G: 24 x 3/4"		F05012434000
G: 25 x 1 1/2"		F05012515000
G: 25 x 1"		F05012510000
G: 25 x 5/8"		F05012506000
G: 25 x 5/8"		F0501250600
G: 26 x 1/2"		F05012605000
G: 26 x 5/8"		F05012606000
G: 27 x 1 1/2"		F05012715000
G: 27 x 1 1/4"		F05012712000
G: 27 x 3/4"		F05012734000
G: 27 x 5/8"		F05012706000
G: 27 x 1/2"		F05012705000
G: 30 x 1/2"		F05013005000

RAMY

