



Blood Administration Set

Al Shifa Medical Co. located at Eastern Region 1st Industrial City, Dammam Kingdom of Saudi Arabia since **1980's**. Al Shifa supplies with a vast range of Medical Devices in the Middle East as well as exports to the different countries of Middle East, Europe, Asia, and Africa.

Al-Shifa supplies a Blood Administration Set, a medical device used for the transfusion of blood or blood components from a blood bag or container to a patient through intravenous infusion. The set is designed to ensure safe, controlled, and sterile delivery of blood products during clinical use.

Al-Shifa Medical Co. supplies different types and configurations of Blood Administration Sets, including variations in tube length, filter type, and optional needle sizes, to meet different clinical requirements.

 **Building a Healthier
Tomorrow**

1st Industrial City | P.O. Box 7917 |
Road 22 | Cross 11 | Dammam-31472 |
Kingdom of Saudi Arabia

RAMY

Product Description

Sterile, Single Use – Blood Administration Set consists of a sharp, non-coring spike for easy insertion into the blood bag port, a transparent drip chamber for visual monitoring of flow, an integrated filter to prevent the passage of clots and debris, and a flexible tube for smooth and controlled blood flow. All materials in contact with blood are non-toxic, pyrogen-free, and suitable for medical use.

Product Intended Use

Sterile, Single Use – Blood Administration Set is intended for the single-use administration of whole blood and blood components from a blood bag to a patient via intravenous infusion under the supervision of a healthcare professional.

Product Features

- Sharp, non-coring spike compatible with blood bags and bottles
- Transparent drip chamber for flow visualization
- Integrated blood filter to trap clots and debris
- Flexible, kink-resistant tubing
- Secure leak-proof connectors
- Available in different tube lengths and configurations
- Suitable for use in adult and pediatric patients
- Single Use, Disposable
- Sterile

Product Types and Variants

Blood Administration Set for Blood Bags

Needle size	Tube length	Catalogue #
W/O needle		FBLDT000010
18G × 1.5"	150 cm	FBAS1815150
18G × 1.5"	180 cm	FBAS1815180
18G × 1.5"	200 cm	FBAS1815200



Blood Collection Set for Blood Bottle with Vented spike

Tube Length	Catalogue #
200 cm	FBLDTB000010

